

Blue Box

Art.N°: KÜ-10044



Advantages

- No expensive cooling logistic necessary.
- Keeping of optimal temperature.
- Reliability for electronic power outage and defects.
- Reusable.
- Stackable.
- Sealable.



3 L cooling elements available as:



- 4°C accu, fluid blue, fluid-cristal-thermometer 0°C–12°C, shows the actual temperature of the element.
For the transport of: red cells (Art.N°: KÜ-10045)
- 22°C accu, fluid green, fluid-cristal-termometer 10°C–40°C, shows the actual temperature of the element
For the transport of: platelets, whole-blood, cord-blood (Art.N°: KÜ-10023)
- - 30°C accu, fluid yellow, for the transport of: fresh frozen plasma

Position of 3 L cooling element:

4 positions above the bloodbags are possible



22 °C



Technical data

Capacity	30 liter
Outer volume (mm)	464 x 368 x 340
Inner volume (mm)	470 x 285 x 230
Material	Expanded

