

Donor lounge Comfort

Available colours

Dialysis/Oncology/Transfusions



Adjustable foot structure

Handset



Comfortable donor lounge for blood banks

- Material: leatherette (different colours), high quality stuffing with perfect density
- Height-adjustable, very comfortable pillow
- Moveable back and seat structure with handset
- Manual or automatic feet section (automation optional)
- Incorporation of 3 motors (4 if feet section is motorized), handset control
- Emergency button (in any place but recommended on the base)
- Safe maximum load of 200 kg
- Removeable armrests, adjustable up to 10° inclination and possibility of lateral and vertical swiveling in order to facilitate access of the patient
- Easy to clean and maintain



Central brake and wheels to move in every direction



Moveable armrest



Shock position

Options

- Additional motor for foot rest adjustment
- Side table (colour black)
- Paper roll holder Inox

Technical data

Dimensions	Seat width 54 cm
	Chair width 88 cm (incl. armrest)
	Total length 206 cm
	Total weight 80 kg
Adjustment	Height 55-75cm
	Shock-position -12°, in 12 seconds
	Pillow adjustable in height
	Armrestadjustment up to 10° inclination
Electrical details	Power Supply: 230v / 50Hz
	On demand: 110v / 50 Hz
	Protection categorie: IP54
	Protection for handset: IP66