

LUXOmatic



- Double press system

High quality system

- Colormatic sensor system for best differentiate between red cells, platelets and plasma
- Long life mechanic
- Ergonomic design



- Sealing heads and clamping for all configurations of Blood bag systems

- Handles for easy lifting of LUXOmatic



Application features

- Flexible system for all types of bags
- Distance diagnose tool integrated- worldwide access for producer and dealer for distance checking and repairing of the machine
- Self adjustable sealing clamps
- Elaborated double press system " precise air system"
- Touch panel display
- 6 sensors for buffy coat control
- 5 valve system for TT bag
- TT 2/3/4 system
- TB 3 /4 system
- Platelets from BC
- Platelets from PRP



- Sensor system
- „Colormatic“ in Multiclamp
- High quality

LUXOmatic

High quality of blood components

- Red Cells – treaded very gentle
- Plasma – cell free because of best sensoric system
- Platelets – recovery rate is very high because of the specific sensoric



- Clamped tube



- Sealing clamp and holding clamp

- Self adjustable sealing clamp for all kinds of blood bags



Easy to handle and operate

- Using a clear organized touch screen
- Price-ergonomic because of two systems in one housing



- Bottom RC bag with permanent LABO balance



- Deposit area to put pencil, clamps, etc. on

Technical data

Power supply	110V-240V / 50/60HZ
Weight	42 kg
Dimensions (cm)	60 x 45 x 50
Power consumption	300 VA or 150 W