

Mobile donor lounge DL 100



The perfect lounge ...



Shock position
Gas pressure supported

Features

- **Portable:** The donor lounge is foldable.
- **Shock position:** A collapsed donor can be set into a shock position supported by a gas tube.
- **Donor slip:** Easy to clean and disinfect. The material has a net structure and is comfortable for the donor.
- **Head rest:** The upholstery of the lounge is integrated in the unit.
- **Arm rest:** The armrest is adjustable on both sides.

... for every blood donation center



Easy to carry
on rolls



Adjustable
armrest



Technical data

| | |
|-------------|--|
| Material: | Stainless steel (round pipe) To prevent scratches and accidents |
| Weight: | 15 kg |
| Dimensions: | fold up: 125 x 55 x 22 cm |