

Spareparts for Sepamatic SL

Name	Code
PC104 Computer board (Wafer)	SM-10007
Disc on chip 256 MB	SM-10019
Power supply 5V,-12V	SM-10044
AD converter board (PC104 AD DIO board)	SM-10030
Step motor board	SM-10031
Interface board (IO)	SM-10029
RF multiplexer board	SM-10032
RF generator	SM-10028
RF cable	SM-10020
RF- Generator + Mux (CSU unit)	SM-10004
Weldinghouse (sealing head)	SM-10005
Cylinder valves set	SM-10045
Double pressure regulator	SM-10142
Valve for double pressure regulator	SM-10143
Pressure sensing switch	SM-10144
Manometer	SM-10093
Close sensor mechanism (End volume switch)	SM-10046
Hall sensor kit	SM-10047
Sensor CX (M-sensors)	SM-10006
Autoteaching sensors (set complete)	SM-10094
Press A pressure plate	SM-10056
Press B moving pressure plate	SM-10133
Press B stable pressure plate	SM-10057
Power transformer	SM-10145

Fuse 6.3A (10 pieces)	SM-90110
LOAD CELL PW2AC3	BA-10047
Bar code reader	SM-10015
LCD panel with keyboard	SM-10014
Keyboard emulator board	SM-10033
Cooler 24 V	SM-10123
Cooler 5V	SM-10122
Pin for balance 1	SM-10082
Pin for balance 2	SM-10081
Pin for balances 3 and 4	SM-10083
Pneumatic cylinder Press A	SM-10146
Pneumatic cylinder Press B	SM-10147
Balance amplifier board	SM-10034
Programming PC-cable SiColMux unit	SM-10011
Stepmotor	SM-10008
Shielded window for LCD	SM-10016
LCD Display	SM-10012
WLAN- Card	SM-90043
Compressor	SM-90019
Flowbreak (Type 1 – Macopress) version 2006	SM-10002
Flowbreak (Type 2 – for clamp A)	SM-10003
PSV2 Board	SM-10036
Power module with transformator (System), lower main supply	SM-10039
Cap for clamp A	SM-10110
Cap for clamp B	SM-10112
Cap for clamp C	SM-10114
Cap for clamp D	SM-10116
Cap for clamp E	SM-10118
Cap for clamp F	SM-10120

Ring for clamp A	SM-10111
Ring for clamp B	SM-10113
Ring for clamp C	SM-10115
Ring for clamp D	SM-10117
Ring for clamp E	SM-10119
Ring for clamp F	SM-10121
Rubber small (for clamps)	SM-10095
Rubber large (for clamps)	SM-10096
Clamp holder fix	SM-10091
Clamp rotating holder movalbe	SM-10092
Sealing ball fixation (bolt)	SM-10090
Breakaway	SM-10059
Holder barcode reader	SM-10038
Peel for SAGM 35 cm	SM-10085
Peel for SAGM 45 cm	SM-10086
Peel for SAGM 50 cm	SM-10099
Sealing electrode full line	SM-10080
Standard packing	SM-10141
Black transport box	SM-90039