

Sterile DMSO

Cryopreservation solutions



Sterile, filtered DMSO is available in a variety of packages and concentrations, and in volumes designed for lab use. 100 % and buffered DMSO solutions are packaged in polymeric vials and syringes; both packages are water-clear, and shatterresistant. Lightweight polymer vials and syringes are free of ions and heavy metals, with a shelf life meeting GLP requirements.




Vials with white caps contain 100 % DMSO and are available in 5 ml, 10 ml, 50 ml, 70 ml and 100 ml volumes with 'flip/tear up' caps

Green top vials have 7.2 ml of 55 % DMSO/ 5 % Dextran-40, also with 'flip/tearup' caps.



Sterile docking syringes

DMSO filled, sterilized syringes are fitted with a stopcock and are shipped disassembled for safety. To use, merely screw the plunger in and connect to your system. They are available in different configurations for easy, aseptic transfer: a length of SCD tubing (shown), a DMSOresistant bag spike or a stopcock cap only.



Sterile DMSO

Product information

- Vial, 5 ml, 100% DMSO, white 'flip/tear up' cap
- Vial, 10 ml, 100% DMSO, white 'flip/tear up' cap
- Vial, 50 ml, 100% DMSO, white 'flip/tear up' cap
- Vial, 70 ml, 100% DMSO, white 'flip/tear up' cap
- Vial, 100 ml, 100% DMSO, white 'flip/tear up' cap
- Syringe 10 ml, 100% DMSO, with stopcock and 4 mm SCD tubing
- Syringe 10 ml, 100% DMSO, with stopcock and cover only
- Syringe 10 ml, 100% DMSO, with stopcock and DMSO resistant bag spike
- Syringe 50 ml, 100% DMSO, with stopcock and 4 mm SCD tubing
- Syringe 50 ml, 100% DMSO, with stopcock and cover only
- Syringe 50 ml, 100% DMSO, with stopcock and DMSO resistant bag spike



55 % DMSO, Buffered with 5 % Dextran-40

- Vial, 7,2 ml, 55 % DMSO and 5 % Dextran-40, green 'flip/tear up' cap
- Syringe, 10 ml, 55 % DMSO and 5 % Dextran-40, with stopcock and 4 mm SCD tubingp
- Syringe, 10 ml, 55 % DMSO and 5 % Dextran-40, with stopcock and cover onlyp
- Syringe, 10 ml, 55 % DMSO and 5 % Dextran-40, with stopcock and DMSO resistant bag spike

- Content: 100 % DMSO Ph.Eur., Or 55% DMSO with 5 % w/v Dextran-40
- Storage: 100 % DMSO - 18 to 25 °C; DMSO/Dextran - 2 to 8°C
- Packaging: Vials and syringes - 12 to a case